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## GB | Support Belt - Work, incl. 4 tool loops

### 1 | Delivery

1 x Support Belt - Work, incl. 4 tool loops

Read the instructions for use completely and carefully!

### 2 | Intended purpose

By external pressure (compression), the abdominal and back muscles should be activated and supported. A stabilization and relief of the lumbar spine and the surrounding musculature can have a relaxing effect and therefore a pain-relieving effect.

Supports the spine and promotes an upright posture.

- Disc problems
- Sciatica
- Lumbago
- Pain from overstrain

For private use

### 3 | Contraindications / Not suitable

In the case of the following medical conditions, it is only recommended to wear this article after consultation with your doctor:

- Skin problems/injuries to the part of the body where the support is worn, especially if there are signs of inflammation, swelling, redness and over-heating
- allergic reactions to the material
- Sensation/movement disorders e.g., associated with diabetes
- open wounds and internal injuries
- allergic reactions to the material
- Lymphatic runoff disorders or unclear soft tissue swelling even outside the scope of application.
- For small children or persons who can not be noticeable when the product sits uncomfortably

### 4 | Notes on safe use of the product

Side effects affecting the whole organism are not yet known. Proper application is required. In case of improper use, product liability is excluded.

- If the bandage is uncomfortable, remove it immediately!
- Wear bandage only by the hour, do not wear permanently!
- If your symptoms worsen or you notice any unusual changes, please contact your doctor immediately.

**Do not use:** while sleeping

### 5 | Application

**Applying support belt:** Fig. 1: First loosen the hook and loop fasteners of the two extra pull straps from the support belt. Hold the support belt with the hook and loop closure (inside) in the right hand and guide it around the back. Pull the two straps sideways away from the body. Fig. 2: Then close the hook and loop fastener at the front at any point on the bandage. Fig. 3: To attach the two extra strings proceed the same way. Fig. 4+5: Cross the two suspenders at the back and attach to the front of the support belt. Fig. 6: Attach the tool loops to the support belt with hook and loop closures

### 6 | Cleaning

Cleaning is done by household laundry, according to care instructions. Close hook and loop fasteners during washing to avoid damaging other items of laundry.

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